

SYNTHETIC POROUS CALCIUM PHOSPHATE GRANULES FOR BONE SUBSTITUTES

20th EAO Annual Scientific Congress, October, 12-15, Athens 2011

Calcium phosphates (CaP) are suitable as bone substitutes due to their biocompatible, bioactive and osteoconductive characteristics. Hydroxyapatite (HA), β -tricalcium phosphate (β TCP) and their combination (HA/ β TCP) are the most used compositions for bone substitutes. In addition, tetracalcium phosphate (TTCP) is a bioceramic that shows in vivo solubility similar to β TCP, but with increased mechanical strength [1]. In the present study are characterized different porous granules suitable as bone lacunae fillers for non load-bearing biomedical applications. Spherical morphology is better than the irregular morphology as the former has greater flexibility in filling irregular lacunae and the latter can be related to inflammatory processes. [2] The combination, HA/TTCP, was attempted and was compared with HA, β TCP and HA/ β TCP manufactured with the same technique. Size control (300-1200 μ m), composition (pure and biphasic calcium phosphates) and internal micro-porosity modulation (0.1-10 μ m) are possible with the presented process. These parameters affect specific surface area and final biomaterial solubility.

MATERIALS & METHODS

CaP materials were supplied in four composition (HA, β TCP, HA/ β TCP, HA/TTCP) in porous granules form (Eurocoating SpA).

Characterization

The morphology, surface roughness, internal porosity and average granules pore size were evaluated by SEM (JSM-5500, Jeol). For the evaluation of the specific surface area (SSA) BET (Nitrogen ASAP 2010 Micromeritics) and Hg-porosimetry analyses (CE 2000) were performed. The density was calculated with He-pycnometry (1035 Micromeritics).

XRD patterns were recorded on Rigaku Dmax III diffractometer by using Cu ka radiation, operating at 40 kV and 30 mA. For FT-IR analysis, granules were reduced to powder and analyzed as KBr pellets with Thermo Nicolet Avatar 330. The solubility of each type of granule was evaluated through dissolution rate and dissolution behavior tests by measuring calcium and phosphate release in TRIS solution at pH 7.3.

The cytotoxicity test was monitored toxic effects of granules on a cell line of mice fibroblasts Balb/c 3T3 in the cell culture under EN ISO 10933-5.

RESULTS & DISCUSSION

Morphology and Microporosity

The round morphology and an example of sub-millimeter dimension for each type of CaP granules are shown in Figure 1.




Fig. 1 - SEM micrographs of spherical a) HA, b) β TCP, c) HA/ β TCP ad d) HA/TTCP granules




Fig. 2 - SEM micrographs of a) HA, b) β TCP, c) HA/ β TCP and d) HA/TTCP granules surface

SEM pictures shown in Figure 2 allow to point out the surface roughness. HA and HA/ β TCP granules have similar homogenous superficial morphology, while β TCP and HA/TTCP spheres are less smooth. It is useful to recall that surface roughness is important for improving idrophilicity and initial adsorption of physiological fluids and blood. Figure 3 shows a comparison of CaP granules section that allows to analyze the internal microporosity.





Fig. 03 - SEM micrographs of a) HA, b) β TCP, c) HA/ β TCP and d) HA/TTCP granules section.

For HA and HA/ β TCP granules micropores are homogeneously distributed and very small (0,1-1 μ m). β TCP and HA/TTCP spheres show bigger pore (1-10 μ m), but the distribution is less regular.

Mineralogical composition

XRD analysis (Fig. 4 and Tab. 1) revealed the final composition of CaP granules; the sum of extraneous phases (CaO for HA, TTCP for HA/ β TCP, α TCP and C2P2O7 for β TCP granules) was less than 5 wt%, thus fulfilling the requirements of ISO 13779-3 norm.



Wwt	HA	β TCP	TTCP	CaO	Ca ₂ P ₂ O ₇	α TCP
β TCP	0,0	99,0	0,0	0,0	0,0	0,0
HA/TTCP	99,8	0,0	49,4	0,0	0,0	0,0
HA/ β TCP	49	49,8	1,1	0,0	0,0	0,0
HA	99,3	0,0	0,0	0,2	0,0	0,0

Figure 4. XRD pattern of SEM micrographs of a) HA, b) HA/ β TCP, c) HA/TTCP and d) β TCP granules

Table 1. Quantitative analysis of CaP granules and detection of foreign phases

FT-IR spectra of crushed CaP granules confirmed the presence of typical bands of calcium phosphate (PO₄³⁻) and the absence of CO₃²⁻, this indicating that all organic phases are burnt out during sintering (Fig. 5 and Fig. 6).




Figure 5. FT-IR spectra of a) HA and b) HA/TTCP granules.




Figure 6. FT-IR spectra of a) HA/ β TCP and b) β TCP granules.

BET , Hg-porosimetry and He-picnometry analyses

The highest SSA was obtained for HA and HA/ β TCP granules that have more small micropores than β TCP. HA/TTCP spheres are the more dense and have the lowest SSA. Table (2)

Solubility test

Ca²⁺ (Fig. 7) and PO₄³⁻ concentration released from granules in TRIS solution referred to 60 days confirmed that HA granules are the less resorbable even with the highest SSA. β TCP granules were also found slowly resorbable due to low SSA. Granules with biphasic compositions are dissolved faster than monophasic ones. Moreover, between biphasics, HA/ β TCP (higher SSA) showed to be more resorbable than HA/TTCP. Finally biphasic HA/TTCP granules resulted more resorbable than monophasic granules in spite of less favorable physical characteristics.

Sample	SSA, BET (m ² /g)	SSA, Hg porosimetry (m ² /g)	Average pore size (μ m)	Porosity (%)	Apparent density (g/cm ³)
HA	3,92	3,70	0,19	47,59	0,64
β TCP	0,67	0,93	3,40	55,76	0,51
HA/ β TCP	2,51	2,36	0,29	49,62	0,47
HA/TTCP	0,48	0,30	0,25	11,11	0,88

Table 2. SSA , average pore size, % porosity and density data for all families.




Figure 6. FT-IR spectra of a) HA/ β TCP and b) β TCP granule

Cytotoxicity test

Cells showed a vitality trend confirming non cytotoxic all the granules composition (i.e. HA; β TCP ; HA/ β TCP; HA/TTCP).

CONCLUSIONS

In this work CaP granules with controlled size, microporosity and different chemical composition were obtained by droplet extrusion method. All granules are spherical with porous surface and internal microstructure can be modulated to favor or limit the dissolution of bone substitute. In vivo test are in progress.

ACKNOWLEDGEMENTS

This work is part of “**CaP Project**” co-sponsored by **Provincia Autonoma di Trento**.

REFERENCES

- [1] C.P.A.T Klein, Studies of the solubility of different calcium phosphate ceramic particles in vitro, *Biomaterials* 11 (1990) p.509
- [2] C. C. Ribeiro · C. C. Barrias · M. A. Barbosa, Preparation and characterization of calcium-phosphate porous microspheres with a uniform size for biomedical applications, *J Mater Sci: Mater Med* (2006) 17: 455–463